

SINGLE-USE FACE SHIELD

- Optically clear, 7 mil thick film
- Full size: Shield length 9"
- Includes foam headband and elastic webbing for comfortable, yet secure fit
- 10 face shields/pack
- USA Made



Dynatomy Single-Use Face Shields are intended to protect personnel against splashing liquids/sprays.

INSTRUCTIONS FOR USE:

- 1. Remove face shield from package.
- 2. Place foam front against your forehead.
- 3. Secure elastic band around the back of your head.
- 4. Discard after single use.
- Please handle your face shield carefully to prevent scratching the shield lens. Do not place or lay the face shield directly on a flat surface as it may scratch the lens.



WE'D LIKE TO HEAR YOUR FEEDBACK.

Please complete a user survey at: daddario.com/shield-survey

Please share your photos on Instagram with the hashtags #daddario and #wewillplayon



Warning: Dynatomy Single-Use Face Shields are lightweight for use in certain light-duty applications, such as medical and dental use. Dynatomy Single-Use Face Shields are not designed to be impact resistant and do not replace ANSI Z87 protective devices, including safety goggles or goggles when needed. Dynatomy Single-Use Face Shields act as a lightweight barrier providing certain secondary protection to help protect against body fluid, sprays or splatter.



Single-Use







DFSM-1-10





D'Addario & Co. 595 Smith Street Farmingdale, NY 11735, USA



06/01/2020

Dynatomy and D'Addario are registered and/ or common law trademarks of D'Addario & Company, Inc., or its affiliates in the United States and/or other countries. Copyright 2020. All rights reserved.

